

产品名称: FABP5 Rabbit Monoclonal Antibody
产品货号: AMRe87021



产品概述 (Summary)

产品名称 (Production Name)	FABP5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FABP5
别名 (Alternative Names)	EFABP; KFABP; E-FABP; PAFABP; PA-FABP
基因 ID (Gene ID)	2171
蛋白 ID (SwissProt ID)	Q01469.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:500-1:2000,ICC/IF 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:15 kDa; Observed MW:15 kDa

研究背景 (Background)


产品名称: FABP5 Rabbit Monoclonal Antibody
产品货号: AMRe87021



This gene encodes the fatty acid binding protein found in epidermal cells, and was first identified as being upregulated in psoriasis tissue. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. FABPs may play roles in fatty acid uptake, transport, and metabolism. Polymorphisms in this gene are associated with type 2 diabetes. The human genome contains many pseudogenes similar to this locus.[provided by RefSeq, Feb 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, T24 cells using AMRe87021 at 1:1000.

注意事项 (Note)

For research use only .