

产品名称: EVPL Rabbit Monoclonal Antibody
产品货号: AMRe87023



产品概述 (Summary)

产品名称 (Production Name)	EVPL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EVPL
别名 (Alternative Names)	EVPK
基因 ID (Gene ID)	2125
蛋白 ID (SwissProt ID)	Q92817.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:232 kDa; Observed MW:232 kDa

研究背景 (Background)


产品名称: EVPL Rabbit Monoclonal Antibody
产品货号: AMRe87023



This gene encodes a member of the plakin family of proteins that forms a component of desmosomes and the epidermal cornified envelope. This gene is located in the tylosis oesophageal cancer locus on chromosome 17q25, and its deletion is associated with both familial and sporadic forms of oesophageal squamous cell carcinoma. Patients suffering from the autoimmune mucocutaneous disorder, paraneoplastic pemphigus, develop antibodies against the encoded protein.
[provided by RefSeq, Jul 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, T24 cells using AMRe87023 at 1:1000.

注意事项 (Note)

For research use only .