

产品名称: DARS Rabbit Monoclonal Antibody
产品货号: AMRe87026



产品概述 (Summary)

产品名称 (Production Name)	DARS Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DARS
别名 (Alternative Names)	HBSL; aspRS
基因 ID (Gene ID)	1615
蛋白 ID (SwissProt ID)	P14868.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:57 kDa; Observed MW:57 kDa

研究背景 (Background)


产品名称: DARS Rabbit Monoclonal Antibody
产品货号: AMRe87026



This gene encodes a member of a multienzyme complex that functions in mediating the attachment of amino acids to their cognate tRNAs. The encoded protein ligates L-aspartate to tRNA(Asp). Mutations in this gene have been found in patients showing hypomyelination with brainstem and spinal cord involvement and leg spasticity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Jurkat, K562 cells using AMRe87026 at 1:1000.

注意事项 (Note)

For research use only .