

产品名称: C1orf93 Rabbit Monoclonal Antibody
产品货号: AMRe87032



产品概述 (Summary)

产品名称 (Production Name)	C1orf93 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C1orf93
别名 (Alternative Names)	C1orf93
基因 ID (Gene ID)	127281
蛋白 ID (SwissProt ID)	Q8TBF2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:21 kDa; Observed MW:21 kDa

研究背景 (Background)


产品名称: C1orf93 Rabbit Monoclonal Antibody
产品货号: AMRe87032




Catalyzes the reduction of prostaglandin-ethanolamide H2 (prostamide H2) to prostamide F(2alpha) with NADPH as proton donor. Also able to reduce prostaglandin H2 to prostaglandin F(2alpha) (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from MCF7 cells using AMRe87032 at 1:1000.



Western blot analysis of extracts from Rat kidney, Mouse brain tissue using AMRe87032 at 1:1000.

注意事项 (Note)

For research use only .