

产品名称: INTS10 Rabbit Monoclonal Antibody
产品货号: AMRe87035



产品概述 (Summary)

产品名称 (Production Name)	INTS10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	INTS10
别名 (Alternative Names)	INT10; C8orf35
基因 ID (Gene ID)	55174
蛋白 ID (SwissProt ID)	Q9NVR2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:82 kDa; Observed MW:82 kDa

研究背景 (Background)


产品名称: INTS10 Rabbit Monoclonal Antibody
产品货号: AMRe87035



INTS10 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005 [PubMed 16239144]).[supplied by OMIM, Mar 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, Jurkat cells using AMRe87035 at 1:1000.

注意事项 (Note)

For research use only .