

产品名称: TNFAIP8 Rabbit Monoclonal Antibody
产品货号: AMRe87040



产品概述 (Summary)

产品名称 (Production Name)	TNFAIP8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNFAIP8
别名 (Alternative Names)	NDED; GG2-1; SCCS2; SCC-S2; MDC-3.13
基因 ID (Gene ID)	25816
蛋白 ID (SwissProt ID)	O95379.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:23 kDa; Observed MW:21-23 kDa

研究背景 (Background)


产品名称: TNFAIP8 Rabbit Monoclonal Antibody
产品货号: AMRe87040




Acts as a negative mediator of apoptosis and may play a role in tumor progression. Suppresses the TNF-mediated apoptosis by inhibiting caspase-8 activity but not the processing of procaspase-8, subsequently resulting in inhibition of BID cleavage and caspase-3 activation.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of TNFAIP8 in CEM, K562 cell lysates using TNFAIP8 antibody (1:1000 diluted).



Western blot analysis of extracts from A549 cells using AMRe87040 at 1:1000.

注意事项 (Note)

For research use only .