

产品名称: TRIM24 Rabbit Monoclonal Antibody
产品货号: AMRe87041



产品概述 (Summary)

产品名称 (Production Name)	TRIM24 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRIM24
别名 (Alternative Names)	PTC6; TF1A; TIF1; RNF82; TIF1A; hTIF1; TIF1ALPHA
基因 ID (Gene ID)	8805
蛋白 ID (SwissProt ID)	O15164.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:117 kDa; Observed MW:117 kDa

研究背景 (Background)


产品名称: TRIM24 Rabbit Monoclonal Antibody
产品货号: AMRe87041




The protein encoded by this gene mediates transcriptional control by interaction with the activation function 2 (AF2) region of several nuclear receptors, including the estrogen, retinoic acid, and vitamin D3 receptors. The protein localizes to nuclear bodies and is thought to associate with chromatin and heterochromatin-associated factors. The protein is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains - a RING, a B-box type 1 and a B-box type 2 - and a coiled-coil region. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of TRIM24 in HeLa, A549 cell lysates using TRIM24 antibody (1:1000 diluted).



Western blot analysis of extracts from HeLa cells using AMRe87041 at 1:1000.

注意事项 (Note)

For research use only.