

产品名称: beta I Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe87052



产品概述 (Summary)

产品名称 (Production Name)	beta I Tubulin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	beta I Tubulin
别名 (Alternative Names)	Tubulin beta-1 chain; TUBB1
基因 ID (Gene ID)	81027
蛋白 ID (SwissProt ID)	Q9H4B7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:200-1:500,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:50 kDa; Observed MW:50 kDa

研究背景 (Background)


产品名称: beta I Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe87052



This gene encodes a member of the beta tubulin protein family. Beta tubulins are one of two core protein families (alpha and beta tubulins) that heterodimerize and assemble to form microtubules. This protein is specifically expressed in platelets and megakaryocytes and may be involved in proplatelet production and platelet release. A mutations in this gene is associated with autosomal dominant macrothrombocytopenia. Two pseudogenes of this gene are found on chromosome Y.[provided by RefSeq, Jul 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from RAW264.7, HeLa cells and Rat brain tissue using AMRe87052 at 1:1000.

注意事项 (Note)

For research use only .