

产品名称: DEC-205 Rabbit Monoclonal Antibody
产品货号: AMRe87055



产品概述 (Summary)

产品名称 (Production Name)	DEC-205 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DEC-205
别名 (Alternative Names)	LY75; CD205; Ly-75; CLEC13B; gp200-MR6
基因 ID (Gene ID)	100526664
蛋白 ID (SwissProt ID)	O60449.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:198 kDa; Observed MW:198 kDa

研究背景 (Background)


产品名称: DEC-205 Rabbit Monoclonal Antibody
产品货号: AMRe87055



This locus represents naturally occurring read-through transcription between the neighboring lymphocyte antigen 75 (LY75) and CD302 molecule (CD302) genes. Alternative splicing results in multiple transcript variants encoding fusion products that share sequence identity with each individual gene product. [provided by RefSeq, Nov 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse thymus tissue using AMRe87055 at 1:1000.

注意事项 (Note)

For research use only .