

产品名称: Clusterin Rabbit Monoclonal Antibody
产品货号: AMRe87059



产品概述 (Summary)

产品名称 (Production Name)	Clusterin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Clusterin
别名 (Alternative Names)	CLI; AAG4; APOJ; CLU1; CLU2; KUB1; SGP2; APO-J; SGP-2; SP-40; TRPM2; TRPM-2; NA1/NA2
基因 ID (Gene ID)	1191
蛋白 ID (SwissProt ID)	P10909.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:53 kDa; Observed MW:35,65 kDa

产品名称: Clusterin Rabbit Monoclonal Antibody
产品货号: AMRe87059




研究背景 (Background)

The protein encoded by this gene is a secreted chaperone that can under some stress conditions also be found in the cell cytosol. It has been suggested to be involved in several basic biological events such as cell death, tumor progression, and neurodegenerative disorders. Alternate splicing results in both coding and non-coding variants.[provided by RefSeq, May 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human fetal lung tissue using AMRe87059 at 1:1000.

注意事项 (Note)

For research use only .