

产品名称: TRAX Rabbit Monoclonal Antibody
产品货号: AMRe87072



产品概述 (Summary)

产品名称 (Production Name)	TRAX Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAX
别名 (Alternative Names)	TRAX
基因 ID (Gene ID)	7257
蛋白 ID (SwissProt ID)	Q99598.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:10000, ICC/IF 1:100-1:200, IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:33 kDa; Observed MW:33 kDa

研究背景 (Background)


产品名称: TRAX Rabbit Monoclonal Antibody
产品货号: AMRe87072




This gene encodes a protein which specifically interacts with translin, a DNA-binding protein that binds consensus sequences at breakpoint junctions of chromosomal translocations. The encoded protein contains bipartite nuclear targeting sequences that may provide nuclear transport for translin, which lacks any nuclear targeting motifs. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of TRAX in K562, Rat Brain, C6, HeLa cell lysates using TRAX antibody (1:1000 diluted).



Western blot analysis of extracts from A431 cells and Rat heart tissue using AMRe87072 at 1:2000.

注意事项 (Note)

For research use only.