

产品名称: SNF2L Rabbit Monoclonal Antibody
产品货号: AMRe87073



产品概述 (Summary)

产品名称 (Production Name)	SNF2L Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNF2L
别名 (Alternative Names)	SWI; ISWI; SWI2; SNF2L; SNF2L1; SNF2LB; SNF2LT; hSNF2L; NURF140
基因 ID (Gene ID)	6594
蛋白 ID (SwissProt ID)	P28370.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:123 kDa; Observed MW:123 kDa

研究背景 (Background)


产品名称: SNF2L Rabbit Monoclonal Antibody
产品货号: AMRe87073



This gene encodes a member of the SWI/SNF family of proteins. The encoded protein is an ATPase which is expressed in diverse tissues and contributes to the chromatin remodeling complex that is involved in transcription. The protein may also play a role in DNA damage, growth inhibition and apoptosis of cancer cells. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from SH-SY5Y cells using AMRe87073 at 1:1000.

注意事项 (Note)

For research use only .