

产品名称: Mannan Binding Lectin Rabbit Monoclonal Antibody
产品货号: AMRe87077



产品概述 (Summary)

产品名称 (Production Name)	Mannan Binding Lectin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Mannan Binding Lectin
别名 (Alternative Names)	MBL; MBP; MBP1; MBPD; MBL2D; MBP-C; COLEC1; HSMBPC
基因 ID (Gene ID)	4153
蛋白 ID (SwissProt ID)	P11226.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:26 kDa; Observed MW:26 kDa

产品名称: Mannan Binding Lectin Rabbit Monoclonal Antibody
产品货号: AMRe87077



研究背景 (Background)

This gene encodes the soluble mannose-binding lectin or mannose-binding protein found in serum. The protein encoded belongs to the collectin family and is an important element in the innate immune system. The protein recognizes mannose and N-acetylglucosamine on many microorganisms, and is capable of activating the classical complement pathway. Deficiencies of this gene have been associated with susceptibility to autoimmune and infectious diseases. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse liver tissue using AMRe87077 at 1:1000.

注意事项 (Note)

For research use only .