

产品名称: GBF1 Rabbit Monoclonal Antibody
产品货号: AMRe87081



产品概述 (Summary)

产品名称 (Production Name)	GBF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GBF1
别名 (Alternative Names)	ARF1GEF
基因 ID (Gene ID)	8729
蛋白 ID (SwissProt ID)	Q92538.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:50-1:100,ICC/IF 1:200-1:500,FC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:207 kDa; Observed MW:250 kDa

研究背景 (Background)


产品名称: GBF1 Rabbit Monoclonal Antibody
产品货号: AMRe87081



This gene encodes a member of the Sec7 domain family. The encoded protein is a guanine nucleotide exchange factor that regulates the recruitment of proteins to membranes by mediating GDP to GTP exchange. The encoded protein is localized to the Golgi apparatus and plays a role in vesicular trafficking by activating ADP ribosylation factor 1. The encoded protein has also been identified as an important host factor for viral replication. Multiple transcript variants have been observed for this gene. [provided by RefSeq, Dec 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HCT116 cells and Mouse heart tissue using AMRe87081 at 1:5000.

注意事项 (Note)

For research use only .