

产品名称: IRAK2 Rabbit Monoclonal Antibody
产品货号: AMRe87082



产品概述 (Summary)

产品名称 (Production Name)	IRAK2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRAK2
别名 (Alternative Names)	IRAK-2
基因 ID (Gene ID)	3656
蛋白 ID (SwissProt ID)	O43187.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:100-1:200,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:69 kDa; Observed MW:69 kDa

研究背景 (Background)


产品名称: IRAK2 Rabbit Monoclonal Antibody
产品货号: AMRe87082



IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from A431 cells using db15341 at 1:1000.

注意事项 (Note)

For research use only .