

产品名称: SAMHD1 Rabbit Monoclonal Antibody
产品货号: AMRe87084



产品概述 (Summary)

产品名称 (Production Name)	SAMHD1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SAMHD1
别名 (Alternative Names)	DCIP; CHBL2; HDDC1; MOP-5; SBB188
基因 ID (Gene ID)	25939
蛋白 ID (SwissProt ID)	Q9Y3Z3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200,ICC/IF 1:100-1:200,FC 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:72 kDa; Observed MW:72 kDa

研究背景 (Background)


产品名称: SAMHD1 Rabbit Monoclonal Antibody
产品货号: AMRe87084




This gene may play a role in regulation of the innate immune response. The encoded protein is upregulated in response to viral infection and may be involved in mediation of tumor necrosis factor-alpha proinflammatory responses. Mutations in this gene have been associated with Aicardi-Goutieres syndrome. [provided by RefSeq, Mar 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, MCF-7, K562 cells using AMRe87084 at 1:1000.



Western blot analysis of extracts from HepG2 cells using AMRe87084 at 1:5000.

注意事项 (Note)

For research use only .