

产品名称: PGK2 Rabbit Monoclonal Antibody
产品货号: AMRe87088



产品概述 (Summary)

产品名称 (Production Name)	PGK2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PGK2
别名 (Alternative Names)	PGKB; PGKPS; HEL-S-272; dJ417L20.2
基因 ID (Gene ID)	5232
蛋白 ID (SwissProt ID)	P07205.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:45 kDa; Observed MW:45 kDa

研究背景 (Background)


产品名称: PGK2 Rabbit Monoclonal Antibody
产品货号: AMRe87088



This gene is intronless, arose via retrotransposition of the phosphoglycerate kinase 1 gene, and is expressed specifically in the testis. Initially assumed to be a pseudogene, the encoded protein is actually a functional phosphoglycerate kinase that catalyzes the reversible conversion of 1,3-bisphosphoglycerate to 3-phosphoglycerate, during the Embden-Meyerhof-Parnas pathway of glycolysis, in the later stages of spermatogenesis.[provided by RefSeq, May 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse testis tissue using AMRe87088 at 1:1000.

注意事项 (Note)

For research use only .