

产品名称: Phospho-Integrin alpha 4 (Ser1027)Rabbit
Monoclonal Antibody
产品货号: AMRe87112



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Integrin alpha 4 (Ser1027)Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Integrin alpha 4
别名 (Alternative Names)	IA4; CD49D
基因 ID (Gene ID)	3676
蛋白 ID (SwissProt ID)	P13612.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:115 kDa; Observed MW:150,70 kDa

**产品名称: Phospho-Integrin alpha 4 (Ser1027)Rabbit
Monoclonal Antibody**
产品货号: AMRe87112




研究背景 (Background)

The gene encodes a member of the integrin alpha chain family of proteins. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain that function in cell surface adhesion and signaling. The encoded preproprotein is proteolytically processed to generate light and heavy chains that comprise the alpha 4 subunit. This subunit associates with a beta 1 or beta 7 subunit to form an integrin that may play a role in cell motility and migration. This integrin is a therapeutic target for the treatment of multiple sclerosis, Crohn's disease and inflammatory bowel disease. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse spleen tissue using AMRe87112 at 1:1000.

注意事项 (Note)

For research use only .