

产品名称: Integrin beta 3 Rabbit Monoclonal Antibody
产品货号: AMRe87131



产品概述 (Summary)

产品名称 (Production Name)	Integrin beta 3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Integrin beta 3
别名 (Alternative Names)	GT; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16
基因 ID (Gene ID)	3690
蛋白 ID (SwissProt ID)	P05106.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:87 kDa; Observed MW:100 kDa

研究背景 (Background)


产品名称: Integrin beta 3 Rabbit Monoclonal Antibody
产品货号: AMRe87131




The ITGB3 protein product is the integrin beta chain beta 3. Integrins are integral cell-surface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. Integrin beta 3 is found along with the alpha IIb chain in platelets. Integrins are known to participate in cell adhesion as well as cell-surface mediated signalling. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Rat heart tissue using AMRe87131 at 1:1000.



Western blot analysis of extracts from HUVEC cells using AMRe87131 at 1:1000.

注意事项 (Note)

For research use only .