

**产品名称: ARFGEF2 Rabbit Monoclonal Antibody**  
**产品货号: AMRe87134**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | ARFGEF2 Rabbit Monoclonal Antibody     |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB                                     |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                          |
|------------------------|--------------------------|
| 基因名 (Gene Name)        | ARFGEF2                  |
| 别名 (Alternative Names) | BIG2; PVNH2; dJ1164I10.1 |
| 基因 ID (Gene ID)        | 10564                    |
| 蛋白 ID (SwissProt ID)   | Q9Y6D5.                  |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000                           |
| 蛋白分子量 (Molecular Weight) | Calculated MW:202 kDa; Observed MW:202 kDa |

## 研究背景 (Background)


**产品名称: ARFGEF2 Rabbit Monoclonal Antibody**  
**产品货号: AMRe87134**



ADP-ribosylation factors (ARFs) play an important role in intracellular vesicular trafficking. The protein encoded by this gene is involved in the activation of ARFs by accelerating replacement of bound GDP with GTP and is involved in Golgi transport. It contains a Sec7 domain, which may be responsible for its guanine-nucleotide exchange activity and also brefeldin A inhibition. [provided by RefSeq, Jul 2008]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from T24 cells using AMRe87134 at 1:3000.

## 注意事项 (Note)

For research use only .