

**产品名称: Estrogen Related Receptor alpha Rabbit
Monoclonal Antibody**
产品货号: AMRe87150



产品概述 (Summary)

产品名称 (Production Name)	Estrogen Related Receptor alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Estrogen Related Receptor alpha
别名 (Alternative Names)	ERR1; ERRa; ESRL1; NR3B1; ERRalpha
基因 ID (Gene ID)	2101
蛋白 ID (SwissProt ID)	P11474.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:46 kDa; Observed MW:50 kDa

**产品名称: Estrogen Related Receptor alpha Rabbit
Monoclonal Antibody**
产品货号: AMRe87150




研究背景 (Background)

The protein encoded by this gene is a nuclear receptor that is most closely related to the estrogen receptor. This protein acts as a site-specific transcription factor and interacts with members of the PGC-1 family of transcription cofactors to regulate the expression of most genes involved in cellular energy production as well as in the process of mitochondrial biogenesis. A processed pseudogene of ESRRA is located on chromosome 13q12.1. [provided by RefSeq, Jun 2019]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87150 at 1:1000.

注意事项 (Note)

For research use only .