

产品名称: CKS2 Rabbit Monoclonal Antibody
产品货号: AMRe87152



产品概述 (Summary)

产品名称 (Production Name)	CKS2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CKS2
别名 (Alternative Names)	CKSHS2
基因 ID (Gene ID)	1164
蛋白 ID (SwissProt ID)	P33552.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:10 kDa; Observed MW:10 kDa

研究背景 (Background)


产品名称: CKS2 Rabbit Monoclonal Antibody
产品货号: AMRe87152




CKS2 protein binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function. The CKS2 mRNA is found to be expressed in different patterns through the cell cycle in HeLa cells, which reflects specialized role for the encoded protein. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Caco-2,293T cells using AMRe87152 at 1:1000.



Western blot analysis of extracts from T24, C2C12 cells and Mouse small intestine tissue using AMRe87152 at 1:1000.

注意事项 (Note)

For research use only .