

产品名称: SAMD9 Rabbit Monoclonal Antibody
产品货号: AMRe87153



产品概述 (Summary)

产品名称 (Production Name)	SAMD9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SAMD9
别名 (Alternative Names)	NFTC; OEF1; OEF2; DRIF1; C7orf5; MIRAGE
基因 ID (Gene ID)	54809
蛋白 ID (SwissProt ID)	Q5K651.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:50-1:100, FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:184 kDa; Observed MW:184 kDa

研究背景 (Background)


产品名称: SAMD9 Rabbit Monoclonal Antibody
产品货号: AMRe87153



This gene encodes a sterile alpha motif domain-containing protein. The encoded protein localizes to the cytoplasm and may play a role in regulating cell proliferation and apoptosis. Mutations in this gene are the cause of normophosphatemic familial tumoral calcinosis. Alternate splicing results in multiple transcript variants that encode the same protein.[provided by RefSeq, Jul 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of SAMD9 in A549 cell lysates using SAMD9 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .