

产品名称: Adrenodoxin Rabbit Monoclonal Antibody
产品货号: AMRe87155



产品概述 (Summary)

产品名称 (Production Name)	Adrenodoxin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Adrenodoxin
别名 (Alternative Names)	ADX; FDX; LOH11CR1D
基因 ID (Gene ID)	2230
蛋白 ID (SwissProt ID)	P10109.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:19 kDa; Observed MW:14 kDa

研究背景 (Background)


产品名称: Adrenodoxin Rabbit Monoclonal Antibody
产品货号: AMRe87155




This gene encodes a small iron-sulfur protein that transfers electrons from NADPH through ferredoxin reductase to mitochondrial cytochrome P450, involved in steroid, vitamin D, and bile acid metabolism. Pseudogenes of this functional gene are found on chromosomes 20 and 21. [provided by RefSeq, Aug 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Adrenodoxin in 3T3 cell lysates using Adrenodoxin antibody(1:1000 diluted).



Western blot analysis of extracts from Mouse kidney tissue using AMRe87155 at 1:1000.

注意事项 (Note)

For research use only .