

产品名称: SULT2A1/ST2 Rabbit Monoclonal Antibody
产品货号: AMRe87157



产品概述 (Summary)

产品名称 (Production Name)	SULT2A1/ST2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SULT2A1/ST2
别名 (Alternative Names)	HST; ST2; STD; hSTa; DHEAS; ST2A1; ST2A3; DHEA-ST
基因 ID (Gene ID)	6822
蛋白 ID (SwissProt ID)	Q06520.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:500-1:2000,ICC/IF 1:20-1:50,FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:34 kDa; Observed MW:34 kDa

研究背景 (Background)


产品名称: SULT2A1/ST2 Rabbit Monoclonal Antibody
产品货号: AMRe87157



This gene encodes a member of the sulfotransferase family. Sulfotransferases aid in the metabolism of drugs and endogenous compounds by converting these substances into more hydrophilic water-soluble sulfate conjugates that can be easily excreted. This protein catalyzes the sulfation of steroids and bile acids in the liver and adrenal glands, and may have a role in the inherited adrenal androgen excess in women with polycystic ovary syndrome. [provided by RefSeq, Mar 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human fetal liver tissue using AMRe87157 at 1:1000.

注意事项 (Note)

For research use only .