

产品名称: DAPK2 Rabbit Monoclonal Antibody
产品货号: AMRe87173



产品概述 (Summary)

产品名称 (Production Name)	DAPK2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DAPK2
别名 (Alternative Names)	DRP1; DRP-1
基因 ID (Gene ID)	23604
蛋白 ID (SwissProt ID)	Q9UIK4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,FC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:43 kDa; Observed MW:43 kDa

研究背景 (Background)


产品名称: DAPK2 Rabbit Monoclonal Antibody
产品货号: AMRe87173



This gene encodes a protein that belongs to the serine/threonine protein kinase family. This protein contains a N-terminal protein kinase domain followed by a conserved calmodulin-binding domain with significant similarity to that of death-associated protein kinase 1 (DAPK1), a positive regulator of programmed cell death. Overexpression of this gene was shown to induce cell apoptosis. It uses multiple polyadenylation sites. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of DAPK2 in Rat spleen cell lysates using DAPK2 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .