

产品名称: IQGAP2 Rabbit Monoclonal Antibody
产品货号: AMRe87181



产品概述 (Summary)

产品名称 (Production Name)	IQGAP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IQGAP2
别名 (Alternative Names)	Ras GTPase-activating-like protein IQGAP2
基因 ID (Gene ID)	10788
蛋白 ID (SwissProt ID)	Q13576.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:181 kDa; Observed MW:181 kDa

研究背景 (Background)


产品名称: IQGAP2 Rabbit Monoclonal Antibody
产品货号: AMRe87181



This gene encodes a member of the IQGAP family. The encoded protein contains three IQ domains, one calponin homology domain, one Ras-GAP domain and one WW domain. This protein interacts with components of the cytoskeleton, with cell adhesion molecules, and with several signaling molecules to regulate cell morphology and motility. It also acts as a tumor suppressor and has been found to play a role in regulating innate antiviral responses. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2017]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from 293T cells and Rat stomach, Rat liver tissue using AMRe87181 at 1:1000.

注意事项 (Note)

For research use only .