

产品名称: MST3 Rabbit Monoclonal Antibody
产品货号: AMRe87185



产品概述 (Summary)

产品名称 (Production Name)	MST3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MST3
别名 (Alternative Names)	MST3; STK3; MST3B; STE20; HEL-S-95
基因 ID (Gene ID)	8428
蛋白 ID (SwissProt ID)	Q9Y6E0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:49 kDa; Observed MW:59 kDa

研究背景 (Background)


产品名称: MST3 Rabbit Monoclonal Antibody
产品货号: AMRe87185




This gene encodes a serine/threonine protein kinase that functions upstream of mitogen-activated protein kinase (MAPK) signaling. The encoded protein is cleaved into two chains by caspases; the N-terminal fragment (MST3/N) translocates to the nucleus and promotes programmed cells death. There is a pseudogene for this gene on chromosome X. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87185 at 1:1000.



Western blot analysis of extracts from HCT116 cells using AMRe87185 at 1:2000.

注意事项 (Note)

For research use only .