

产品名称: GPCR GPR26 Rabbit Monoclonal Antibody
产品货号: AMRe87206



产品概述 (Summary)

产品名称 (Production Name)	GPCR GPR26 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPCR GPR26
别名 (Alternative Names)	GPR26
基因 ID (Gene ID)	2849
蛋白 ID (SwissProt ID)	Q8NDV2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:38 kDa; Observed MW:38 kDa

研究背景 (Background)


产品名称: GPCR GPR26 Rabbit Monoclonal Antibody
产品货号: AMRe87206



This gene encodes a G protein-couple receptor protein. G-protein-coupled receptors are a large family of membrane proteins that are involved in cellular responses to environmental stimuli, neurotransmitters, and hormones. The encoded protein may play a role in neurodegenerative diseases. Epigenetic silencing of this gene has been observed in gliomas.
[provided by RefSeq, Sep 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from C6 cells using AMRe87206 at 1:1000.

注意事项 (Note)

For research use only .