

产品名称: Torsin A/DYT1 Rabbit Monoclonal Antibody
产品货号: AMRe87213



产品概述 (Summary)

产品名称 (Production Name)	Torsin A/DYT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Torsin A/DYT1
别名 (Alternative Names)	DQ2; DYT1
基因 ID (Gene ID)	1861
蛋白 ID (SwissProt ID)	O14656.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:38 kDa; Observed MW:35 kDa

研究背景 (Background)


产品名称: Torsin A/DYT1 Rabbit Monoclonal Antibody
产品货号: AMRe87213



The protein encoded by this gene is a member of the AAA family of adenosine triphosphatases (ATPases), is related to the Clp protease/heat shock family and is expressed prominently in the substantia nigra pars compacta. Mutations in this gene result in the autosomal dominant disorder, torsion dystonia 1. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from 3T3 cells using AMRe87213 at 1:1000.

注意事项 (Note)

For research use only .