

产品名称: PITX3/PTX3 Rabbit Monoclonal Antibody
产品货号: AMRe87235



产品概述 (Summary)

产品名称 (Production Name)	PITX3/PTX3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PITX3/PTX3
别名 (Alternative Names)	ASMD; ASOD; PTX3; ASGD1; CTPP4; CTRCT11
基因 ID (Gene ID)	5309
蛋白 ID (SwissProt ID)	O75364.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:32 kDa; Observed MW:32 kDa

研究背景 (Background)


产品名称: PITX3/PTX3 Rabbit Monoclonal Antibody
产品货号: AMRe87235



This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. Members of this family act as transcription factors. This protein is involved in lens formation during eye development. Mutations of this gene have been associated with anterior segment mesenchymal dysgenesis and congenital cataracts. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from 293T,U87-MG cells using AMRe87235 at 1:1000.

注意事项 (Note)

For research use only .