

产品名称: Cathepsin Z Rabbit Monoclonal Antibody
产品货号: AMRe87239



产品概述 (Summary)

产品名称 (Production Name)	Cathepsin Z Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cathepsin Z
别名 (Alternative Names)	CTSX
基因 ID (Gene ID)	1522
蛋白 ID (SwissProt ID)	Q9UBR2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:34 kDa; Observed MW:34 kDa

研究背景 (Background)


产品名称: Cathepsin Z Rabbit Monoclonal Antibody
产品货号: AMRe87239



The protein encoded by this gene is a lysosomal cysteine proteinase and member of the peptidase C1 family. It exhibits both carboxy-monopeptidase and carboxy-dipeptidase activities. The encoded protein has also been known as cathepsin X and cathepsin P. This gene is expressed ubiquitously in cancer cell lines and primary tumors and, like other members of this family, may be involved in tumorigenesis. [provided by RefSeq, Oct 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from T24, HeLa cells using AMRe87239 at 1:2000.

注意事项 (Note)

For research use only .