

产品名称: nNOS (neuronal) Rabbit Monoclonal Antibody
产品货号: AMRe87245



产品概述 (Summary)

产品名称 (Production Name)	nNOS (neuronal) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	nNOS (neuronal)
别名 (Alternative Names)	NOS; bNOS; nNOS; IHPS1; N-NOS; NC-NOS
基因 ID (Gene ID)	4842
蛋白 ID (SwissProt ID)	P29475.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:200, FC 1:200-1:500, IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:161 kDa; Observed MW:161 kDa

研究背景 (Background)


产品名称: nNOS (neuronal) Rabbit Monoclonal Antibody
产品货号: AMRe87245



The protein encoded by this gene belongs to the family of nitric oxide synthases, which synthesize nitric oxide from L-arginine. Nitric oxide is a reactive free radical, which acts as a biologic mediator in several processes, including neurotransmission, and antimicrobial and antitumoral activities. In the brain and peripheral nervous system, nitric oxide displays many properties of a neurotransmitter, and has been implicated in neurotoxicity associated with stroke and neurodegenerative diseases, neural regulation of smooth muscle, including peristalsis, and penile erection. This protein is ubiquitously expressed, with high level of expression in skeletal muscle. Multiple transcript variants that differ in the 5' UTR have been described for this gene but the full-length nature of these transcripts is not known. Additionally, alternatively spliced transcript variants encoding different isoforms (some testis-specific) have been found for this gene.[provided by RefSeq, Feb 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse brain tissue using AMRe87245 at 1:1000.

注意事项 (Note)

For research use only .