

产品名称: ZNF160 Rabbit Monoclonal Antibody
产品货号: AMRe87252



产品概述 (Summary)

产品名称 (Production Name)	ZNF160 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF160
别名 (Alternative Names)	F11; HZF5; KR18; HKr18
基因 ID (Gene ID)	90338
蛋白 ID (SwissProt ID)	Q9HCG1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:10-1:100,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:94 kDa; Observed MW:74 kDa

研究背景 (Background)


产品名称: ZNF160 Rabbit Monoclonal Antibody
产品货号: AMRe87252



The protein encoded by this gene is a Kruppel-related zinc finger protein which is characterized by the presence of an N-terminal repressor domain, the Kruppel-associated box (KRAB). The KRAB domain is a potent repressor of transcription; thus this protein may function in transcription regulation. Multiple transcript variants have been found for this gene.
[provided by RefSeq, Apr 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of ZNF160 in mouse liver cell lysates using ZNF160 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .