

产品名称: Phospho-Nucleolin (Thr84) Rabbit Monoclonal Antibody
产品货号: AMRe87256



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Nucleolin (Thr84) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Nucleolin (Thr84)
别名 (Alternative Names)	C23; Nsr1
基因 ID (Gene ID)	4691
蛋白 ID (SwissProt ID)	P19338.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:77 kDa; Observed MW:100 kDa

产品名称: Phospho-Nucleolin (Thr84) Rabbit Monoclonal Antibody
产品货号: AMRe87256




研究背景 (Background)


Nucleolin (NCL), a eukaryotic nucleolar phosphoprotein, is involved in the synthesis and maturation of ribosomes. It is located mainly in dense fibrillar regions of the nucleolus. Human NCL gene consists of 14 exons with 13 introns and spans approximately 11kb. The intron 11 of the NCL gene encodes a small nucleolar RNA, termed U20. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Phospho-Nucleolin (Thr84) in U251MG, HepG2, Hela cell lysates using Phospho-Nucleolin (Thr84) antibody (1:500 diluted).



Western blot analysis of extracts from U2OS, Jurkat, SCC-9 cells using AMRe87256 at 1:3000.

注意事项 (Note)

For research use only.