

产品名称: CRMP1 Rabbit Monoclonal Antibody
产品货号: AMRe87261



产品概述 (Summary)

产品名称 (Production Name)	CRMP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRMP1
别名 (Alternative Names)	DRP1; DRP-1; CRMP-1; DPYSL1; ULIP-3
基因 ID (Gene ID)	1400
蛋白 ID (SwissProt ID)	Q14194.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:1000-1:5000,ICC/IF 1:200-1:500,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:62,74 kDa; Observed MW:62,74 kDa

产品名称: CRMP1 Rabbit Monoclonal Antibody
产品货号: AMRe87261




研究背景 (Background)

This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from C6 cells using AMRe87261 at 1:1000.

注意事项 (Note)

For research use only .