

产品名称: gamma Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe87262



产品概述 (Summary)

产品名称 (Production Name)	gamma Tubulin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	gamma Tubulin
别名 (Alternative Names)	TUBG; GCP-1; CDCBM4; TUBGCP1
基因 ID (Gene ID)	7283
蛋白 ID (SwissProt ID)	P23258.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:51 kDa; Observed MW:50 kDa

研究背景 (Background)

产品名称: gamma Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe87262



This gene encodes a member of the tubulin superfamily. The encoded protein localizes to the centrosome where it binds to microtubules as part of a complex referred to as the gamma-tubulin ring complex. The protein mediates microtubule nucleation and is required for microtubule formation and progression of the cell cycle. A pseudogene of this gene is found on chromosome 7. [provided by RefSeq, Jan 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87262 at 1:1000.

注意事项 (Note)

For research use only .