

产品名称: Phospho-YAP1 (Ser127) Rabbit Monoclonal Antibody
产品货号: AMRe87275



产品概述 (Summary)

产品名称 (Production Name)	Phospho-YAP1 (Ser127) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-YAP1 (Ser127)
别名 (Alternative Names)	YAP; YKI; COB1; YAP2; YAP65
基因 ID (Gene ID)	10413
蛋白 ID (SwissProt ID)	P46937.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:55 kDa; Observed MW:70-75 kDa

产品名称: Phospho-YAP1 (Ser127) Rabbit Monoclonal Antibody
产品货号: AMRe87275




研究背景 (Background)


This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Phospho-YAP1 (Ser127) in HeLa, HeLa+CalyculinA using Phospho-YAP1 (Ser127) antibody.



Western blot analysis of extracts from T24 cells and Rat spleen, Mouse muscle tissue using AMRe87275 at 1:6000.

注意事项 (Note)

For research use only .