

产品名称: SAP97 Rabbit Monoclonal Antibody
产品货号: AMRe87285



产品概述 (Summary)

产品名称 (Production Name)	SAP97 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SAP97
别名 (Alternative Names)	hdlg; DLGH1; SAP97; SAP-97
基因 ID (Gene ID)	1739
蛋白 ID (SwissProt ID)	Q12959.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:101 kDa; Observed MW:101 kDa

研究背景 (Background)


产品名称: SAP97 Rabbit Monoclonal Antibody
产品货号: AMRe87285




This gene encodes a multi-domain scaffolding protein that is required for normal development. This protein may have a role in septate junction formation, signal transduction, cell proliferation, synaptogenesis and lymphocyte activation. A multitude of transcript variants deriving from alternative splicing and the use of multiple alternate promoter have been observed, including some splice variants that may be specific to brain and other tissues. An upstream uORF may regulate translation at some splice variants of this gene. [provided by RefSeq, Sep 2018]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of SAP97 in HeLa,A549,HL-60 using SAP97 antibody(1:1000 diluted)



Western blot analysis of extracts from A431 cells using AMRe87285 at 1:1000.

注意事项 (Note)

For research use only .