

产品名称: Apolipoprotein CIII Rabbit Monoclonal Antibody
产品货号: AMRe87286



产品概述 (Summary)

产品名称 (Production Name)	Apolipoprotein CIII Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Apolipoprotein CIII
别名 (Alternative Names)	APOCIII
基因 ID (Gene ID)	345
蛋白 ID (SwissProt ID)	P02656.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:11 kDa; Observed MW:11 kDa


研究背景 (Background)

产品名称: Apolipoprotein CIII Rabbit Monoclonal Antibody
产品货号: AMRe87286

This gene encodes a protein component of triglyceride (TG)-rich lipoproteins (TRLs) including very low density lipoproteins (VLDL), high density lipoproteins (HDL) and chylomicrons. The encoded protein plays a role in the metabolism of these TRLs through multiple modes. This protein has been shown to promote the secretion of VLDL1, inhibit lipoprotein lipase enzyme activity, and delay catabolism of TRL remnants. Mutations in this gene are associated with low plasma triglyceride levels and reduced risk of ischemic cardiovascular disease, and hyperalphalipoproteinemia, which is characterized by elevated levels of high density lipoprotein (HDL) and HDL cholesterol in human patients. This gene and other related genes comprise an apolipoprotein gene cluster on chromosome 11. [provided by RefSeq, Sep 2017]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human plasma tissue using AMRe87286 at 1:1000.

注意事项 (Note)

For research use only .