

产品名称: MTAP Rabbit Monoclonal Antibody
产品货号: AMRe87289



产品概述 (Summary)

产品名称 (Production Name)	MTAP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MTAP
别名 (Alternative Names)	BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249
基因 ID (Gene ID)	4507
蛋白 ID (SwissProt ID)	Q13126.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC/IF 1:100-1:200,FC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:31 kDa; Observed MW:31 kDa


研究背景 (Background)

产品名称: MTAP Rabbit Monoclonal Antibody
产品货号: AMRe87289


This gene encodes an enzyme that plays a major role in polyamine metabolism and is important for the salvage of both adenine and methionine. The encoded enzyme is deficient in many cancers because this gene and the tumor suppressor p16 gene are co-deleted. Multiple alternatively spliced transcript variants have been described for this gene, but their full-length natures remain unknown. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of MTAP in A549, MCF-7, U251, HL-60 cell lysates using MTAP antibody (1:1000 diluted).



Western blot analysis of extracts from A431 cells and Mouse brain tissue using AMRe87289 at 1:1000.

注意事项 (Note)

For research use only.