

产品名称: TAF1C Rabbit Monoclonal Antibody
产品货号: AMRe87310



产品概述 (Summary)

产品名称 (Production Name)	TAF1C Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAF1C
别名 (Alternative Names)	SL1; TAFI95; TAFI110; MGC:39976
基因 ID (Gene ID)	9013
蛋白 ID (SwissProt ID)	Q15572.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:95 kDa; Observed MW:95 kDa

研究背景 (Background)


产品名称: TAF1C Rabbit Monoclonal Antibody
产品货号: AMRe87310




Initiation of transcription by RNA polymerase I requires the formation of a complex composed of the TATA-binding protein (TBP) and three TBP-associated factors (TAFs) specific for RNA polymerase I. This complex, known as SL1, binds to the core promoter of ribosomal RNA genes to position the polymerase properly and acts as a channel for regulatory signals. This gene encodes the largest SL1-specific TAF. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of TAF1C in K562 cell lysates using TAF1C antibody(1:1000 diluted).



Western blot analysis of extracts from HeLa, T24 cells using AMRe87310 at 1:1000.

注意事项 (Note)

For research use only .