

产品名称: Laminin Rabbit Monoclonal Antibody
产品货号: AMRe87317



产品概述 (Summary)

产品名称 (Production Name)	Laminin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Laminin
别名 (Alternative Names)	AI562206; AW240805; 1110064A23Rik
基因 ID (Gene ID)	23928
蛋白 ID (SwissProt ID)	Q9R0B6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:172 kDa; Observed MW:172 kDa

研究背景 (Background)


产品名称: Laminin Rabbit Monoclonal Antibody
产品货号: AMRe87317



Binding to cells via a high affinity receptor, laminin is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse kidney tissue using AMRe87317 at 1:1000.

注意事项 (Note)

For research use only .