

产品名称: HDLBP Rabbit Monoclonal Antibody
产品货号: AMRe87335



产品概述 (Summary)

产品名称 (Production Name)	HDLBP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HDLBP
别名 (Alternative Names)	HBP; VGL; PRO2900
基因 ID (Gene ID)	3069
蛋白 ID (SwissProt ID)	Q00341.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:142 kDa; Observed MW:142 kDa

研究背景 (Background)


产品名称: HDLBP Rabbit Monoclonal Antibody
产品货号: AMRe87335




The protein encoded by this gene binds high density lipoprotein (HDL) and may function to regulate excess cholesterol levels in cells. The encoded protein also binds RNA and can induce heterochromatin formation. [provided by RefSeq, Mar 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of HDLBP in HeLa, A549, HL-60, U251, U87-MG cell lysates using HDLBP antibody (1:1000 diluted).



Western blot analysis of extracts from HeLa, C2C12 cells and Mouse muscle, Rat heart tissue using AMRe87335 at 1:5000.

注意事项 (Note)

For research use only.