

产品名称: IGF2BP1 Rabbit Monoclonal Antibody
产品货号: AMRe87339



产品概述 (Summary)

产品名称 (Production Name)	IGF2BP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGF2BP1
别名 (Alternative Names)	IMP1; ZBP1; CRDBP; IMP-1; CRD-BP; VICKZ1
基因 ID (Gene ID)	10642
蛋白 ID (SwissProt ID)	Q9NZI8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:64 kDa; Observed MW:64 kDa

研究背景 (Background)


产品名称: IGF2BP1 Rabbit Monoclonal Antibody
产品货号: AMRe87339




This gene encodes a member of the insulin-like growth factor 2 mRNA-binding protein family. The protein encoded by this gene contains four K homology domains and two RNA recognition motifs. It functions by binding to the mRNAs of certain genes, including insulin-like growth factor 2, beta-actin and beta-transducin repeat-containing protein, and regulating their translation. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of IGF2BP1 in K562,3T3 cell lysates using IGF2BP1 antibody(1:1000 diluted).



Western blot analysis of extracts from T24 cells using AMRe87339 at 1:1000.

注意事项 (Note)

For research use only .