

产品名称: Phospho-p130 (Thr986) Rabbit Monoclonal Antibody
产品货号: AMRe87345



产品概述 (Summary)

产品名称 (Production Name)	Phospho-p130 (Thr986) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-p130
别名 (Alternative Names)	Rb2; P130
基因 ID (Gene ID)	5934
蛋白 ID (SwissProt ID)	Q08999.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:100-1:200,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:128 kDa; Observed MW:130 kDa

产品名称: Phospho-p130 (Thr986) Rabbit Monoclonal Antibody
产品货号: AMRe87345




研究背景 (Background)


Key regulator of entry into cell division. Directly involved in heterochromatin formation by maintaining overall chromatin structure and, in particular, that of constitutive heterochromatin by stabilizing histone methylation. Recruits and targets histone methyltransferases KMT5B and KMT5C, leading to epigenetic transcriptional repression. Controls histone H4 'Lys-20' trimethylation. Probably acts as a transcription repressor by recruiting chromatin-modifying enzymes to promoters. Potent inhibitor of E2F-mediated trans-activation, associates preferentially with E2F5. Binds to cyclins A and E. Binds to and may be involved in the transforming capacity of the adenovirus E1A protein. May act as a tumor suppressor.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Phospho-Rb2 p130 (Thr986) in K562 cell lysates using Phospho-Rb2 p130 (Thr986) antibody(1:1000 diluted).



Western blot analysis of extracts from C2C12 cells and Mouse muscle tissue using AMRe87345 at 1:1000.

注意事项 (Note)

For research use only .

产品名称: Phospho-p130 (Thr986) Rabbit Monoclonal
Antibody
产品货号: AMRe87345

