

产品名称: USF1 Rabbit Monoclonal Antibody
产品货号: AMRe87366



产品概述 (Summary)

产品名称 (Production Name)	USF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USF1
别名 (Alternative Names)	UEF; FCHL; MLTF; FCHL1; MLTFI; HYPLIP1; bHLHb11
基因 ID (Gene ID)	7391
蛋白 ID (SwissProt ID)	P22415.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,FC 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:34 kDa; Observed MW:34 kDa

研究背景 (Background)


产品名称: USF1 Rabbit Monoclonal Antibody
产品货号: AMRe87366




This gene encodes a member of the basic helix-loop-helix leucine zipper family, and can function as a cellular transcription factor. The encoded protein can activate transcription through pyrimidine-rich initiator (Inr) elements and E-box motifs. This gene has been linked to familial combined hyperlipidemia (FCHL). Alternative splicing of this gene results in multiple transcript variants. A related pseudogene has been defined on chromosome 21. [provided by RefSeq, Feb 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of USF1 in HeLa,mouse Brain,Ramos using USF1 antibody(1:1000 diluted)



Western blot analysis of extracts from Rat heart tissue using AMRe87366 at 1:5000.

注意事项 (Note)

For research use only .