

产品名称: Ube2L6 Rabbit Monoclonal Antibody
产品货号: AMRe87368



产品概述 (Summary)

产品名称 (Production Name)	Ube2L6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ube2L6
别名 (Alternative Names)	RIG-B; UBCH8
基因 ID (Gene ID)	9246
蛋白 ID (SwissProt ID)	O14933.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:18 kDa; Observed MW:18 kDa

研究背景 (Background)


产品名称: Ube2L6 Rabbit Monoclonal Antibody
产品货号: AMRe87368



The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes (E1s), ubiquitin-conjugating enzymes (E2s) and ubiquitin-protein ligases (E3s). This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is highly similar in primary structure to the enzyme encoded by the UBE2L3 gene. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Molt-4 cells using AMRe87368 at 1:1000.

注意事项 (Note)

For research use only .