

产品名称: Serpin A5 Rabbit Monoclonal Antibody
产品货号: AMRe87376



产品概述 (Summary)

产品名称 (Production Name)	Serpin A5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Serpin A5
别名 (Alternative Names)	PCI; PAI3; PAI-3; PCI-B; PROCI; PLANH3
基因 ID (Gene ID)	5104
蛋白 ID (SwissProt ID)	P05154.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:1000-1:5000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:46 kDa; Observed MW:46-52 kDa

研究背景 (Background)


产品名称: Serpin A5 Rabbit Monoclonal Antibody
产品货号: AMRe87376



The protein encoded by this gene is a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. This gene is one in a cluster of serpin genes located on the q arm of chromosome 14. This family member is a glycoprotein that can inhibit several serine proteases, including protein C and various plasminogen activators and kallikreins, and it thus plays diverse roles in hemostasis and thrombosis in multiple organs. [provided by RefSeq, Aug 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human urine tissue using AMRe87376 at 1:1000.

注意事项 (Note)

For research use only .