

产品名称: TPST2 Rabbit Monoclonal Antibody
产品货号: AMRe87377



产品概述 (Summary)

产品名称 (Production Name)	TPST2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TPST2
别名 (Alternative Names)	TPST-2; TANGO13B
基因 ID (Gene ID)	8459
蛋白 ID (SwissProt ID)	O60704.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:42 kDa; Observed MW:42 kDa

研究背景 (Background)


产品名称: TPST2 Rabbit Monoclonal Antibody
产品货号: AMRe87377



The protein encoded by this gene catalyzes the O-sulfation of tyrosine residues within acidic regions of proteins. The encoded protein is a type II integral membrane protein found in the Golgi body. Alternative splicing produces multiple transcript variants encoding distinct isoforms. [provided by RefSeq, May 2018]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human fetal liver tissue using AMRe87377 at 1:1000.

注意事项 (Note)

For research use only .